

The Buyer's Guide to

Buffer Formulation Equipment

Acquiring the Equipment Needed for Large-Scale Bioproduction





Table of Contents

Introduction	04
Buffer Formulation Equipment	06
Features to Consider	13
Budget & Vendor Considerations	16
Conclusion	21



Introduction

Buffer formulation equipment is integral to your ability to produce biopharmaceuticals at scale. In this buyer's guide, our buffer formulation experts give actionable insights and guidance for getting started.

Biopharmaceuticals are critical in treating various diseases, from cancer to autoimmune disorders. These medical drugs include vaccines, monoclonal antibodies, cell therapies, and other biological products. They are designed to target specific pathways within the body, offering precise and personalized treatment options.

The need for biopharmaceuticals is growing rapidly. As the global population increases, so does the prevalence of complex and chronic diseases, which drives the demand for more effective and innovative treatments.

Simultaneously, advancements in biotechnology are revolutionizing drug development, enabling the creation of highly targeted therapies. These innovations are pushing the industry to produce cutting-edge drugs more efficiently and scale up bioproduction to meet the growing global need. This is advantageous because large-scale bioproduction allows for:



Cost Efficiency

Large-scale bioproduction reduces per-unit costs through economies of scale, making biopharmaceuticals more affordable.



Enhanced Quality Control

Standardized processes and automation in large-scale production improve consistency and reliability.



Faster Distribution

Global production sites enable quicker delivery of medicines to various regions, which minimizes shipping delays and improves accessibility.

How Does Buffer Management Play a Role in Large-Scale Bioproduction?

Buffer management controls and maintains the composition and concentration of buffer solutions essential for large-scale bioproduction. **This process plays a key role in:**

- ▶ **Maintaining optimal pH:** Buffers stabilize the pH of solutions, which is crucial for creating the ideal environment for biochemical reactions and ensuring consistent product quality.
- ▶ **Stabilizing biomolecules:** Proper buffering prevents the degradation and denaturation of sensitive biomolecules, preserving their functionality throughout the production process.
- ▶ **Enhancing reaction efficiency:** Effective buffer preparation ensures biochemical reactions occur under optimal conditions, which improves reaction rates and overall production efficiency.
- ▶ **Facilitating purification processes:** Buffers support chromatographic and other purification methods by maintaining the conditions needed for effective separation and purification of target biomolecules.

Traditional Buffer Preparation vs. Inline Buffer Formulation

Traditional buffer preparation methods involve producing large quantities of buffer solutions and storing them in extensive tank farms until needed. However, if the stored buffer is not used before it deteriorates, it must be purged, and the tanks cleaned and refilled, leading to inefficiencies and additional maintenance costs.

Another common method is the manual preparation of buffer concentrates. These concentrates, which may come in powder form, need to be dissolved and mixed into larger tanks. This process can expose workers to airborne particles from the powders, requiring personal protective equipment (PPE) and increasing the risk of contamination. The manual preparation process can also be labor-intensive and prone to errors.

Addressing these challenges involves exploring more streamlined and automated buffer management solutions to enhance efficiency and safety.

Inline buffer formulation is an ideal solution. This technique involves buffers that are prepared in real time and on demand as needed during the production process rather than in large batches ahead of time. By integrating the buffer preparation directly into the production line, you allow for more flexible, efficient, and scalable operations.

In this buyer's guide, we offer a comprehensive view of inline buffer formulation equipment that can enhance your process. We also provide key considerations for features, budget, and vendors to help you make an educated decision.



Buffer Formulation **Equipment**

The Basics of Downstream Processing for Inline Buffer Formulation

Downstream processes require large volumes of dilute and often complex buffers. Here is how inline buffer formulation plays a role in each step of the downstream process.

Step 1: Cell Harvesting

In the initial step of cell harvesting, prepared buffers help maintain cell integrity and optimize the separation of cells from the culture medium, reducing the risk of premature lysis and ensuring a clean start to downstream processing.

Step 2: Cell Lysis

Moving into cell lysis, buffers provide the necessary conditions for controlled lysis, protecting the stability of the target biomolecules while minimizing the release of unwanted contaminants.

Step 3: Clarification

During clarification, inline buffer systems adjust pH and ionic strength in real time, facilitating the efficient removal of cell debris and preventing product loss.

Step 4: Chromatography

As the process advances to chromatography, inline buffer formulation becomes even more critical. The process delivers precisely tailored buffers for binding, washing, and elution phases, which are essential for the effective separation and purification of the target product. The ability to adjust these buffers dynamically enhances the consistency and quality of the purified product.

Step 5: Filtration

Inline buffer formulation (IBF™) systems support concentration and diafiltration steps by providing fresh buffers that maintain product stability and facilitate buffer exchange, which is crucial for preparing the product for final formulation or further processing.

Step 6: Precipitation

In the precipitation stage, buffers selectively precipitate either the target protein or unwanted impurities by finely tuning the pH and ionic conditions. This controlled precipitation process maximizes yield and purity.

Step 7: Virus Inactivation/Removal

During virus inactivation or removal, inline buffer formulation ensures the product remains stable while effectively eliminating viral contaminants through low pH treatment or filtration.

Step 8: Polishing

Finally, in the polishing step, which aims to achieve the highest purity levels, buffers are customized to fine-tune purification conditions and enhance product stability, ensuring that the final biopharmaceutical product meets stringent regulatory and quality standards.

Across all these stages, inline buffer formulation provides real-time adaptability, precision, and efficiency.

To enhance the inline buffer formulation process, you need equipment that will:

- ▶ Precisely control buffer composition.
- ▶ Automate real-time adjustments.
- ▶ Integrate seamlessly with downstream processing systems.
- ▶ Ensure consistent, high-quality buffer preparation tailored to the specific needs of each production stage.

Let's explore them now.

Innovative Inline Buffer Formulation Equipment

Asahi Kasei Bioprocess offers the following equipment to improve the inline buffer formulation process:

The MOTIV 5 logo features the word "MOTIV" in a bold, blue, sans-serif font, followed by a large, bold number "5". Below the letters "O", "T", and "I" is a decorative horizontal line of blue dots of varying sizes, creating a sense of motion or a buffer layer.The MOTIV plus logo features the word "MOTIV" in a bold, blue, sans-serif font, followed by the word "plus" in a blue, italicized, sans-serif font. Below the letters "O", "T", and "I" is a decorative horizontal line of blue dots of varying sizes, creating a sense of motion or a buffer layer.The MOTIV 3 logo features the word "MOTIV" in a bold, blue, sans-serif font, followed by a large, bold number "3". Below the letters "O", "T", and "I" is a decorative horizontal line of blue dots of varying sizes, creating a sense of motion or a buffer layer.The MOTIV SU logo features the word "MOTIV" in a bold, blue, sans-serif font, followed by the letters "SU" in a bold, blue, sans-serif font. Below the letters "O", "T", and "I" is a decorative horizontal line of blue dots of varying sizes, creating a sense of motion or a buffer layer.



MOTIV[®] 3

MOTIV[®] 3

The [MOTIV 3](#) is a 3-pump inline buffer formulation system that simplifies the creation of dilute and conditioned buffers using up to three separate streams.

A larger pump delivers the diluent (typically WFI or purified water) into the system, while two smaller pumps supply the buffer concentrate and an acid or base modifier.

There are several product options for the [MOTIV 3](#) that offer different flow rates:

Name	Liters per hour
MOTIV [®] 310	Up to 1,000
MOTIV [®] 312	Up to 1,200
MOTIV [®] 325	Up to 2,500
MOTIV [®] 350	Up to 5,000

MOTIV[®] 5


MOTIV[®] 5

The [MOTIV 5](#) is a 5-pump inline buffer formulation system designed to produce complex buffers using up to five different streams.

The system allows for precise buffer formulation through conductivity, pH, and mass flow control.

MOTIV 5 comes with varying flow rates:

Name	Liters per hour
MOTIV [®] 510	Up to 1,000
MOTIV [®] 512	Up to 1,200
MOTIV [®] 525	Up to 2,500
MOTIV [®] 550	Up to 5,000





MOTIV[®]plus

MOTIV[®]plus

If the standard 3- and 5-pump designs don't meet your production needs, Asahi Kasei Bioprocess can build you a system that will.

The [MOTIV[®]plus](#) offers enhanced customization options and supports non-standard flow rate ranges. The system uses specialized materials such as AL6XN for corrosion resistance and accommodates additional pumps for processes that require more than five primary streams.

MOTIVplus can also integrate additional inline monitors, such as refractive index sensors, to provide greater control and precision in buffer formulation.

MOTIV[®] SU

MOTIV[®] SU

The [MOTIV SU](#) is a single-use inline buffer formulation system designed to produce complex buffers on demand with efficiency and precision. It creates buffers all from one pump head without requiring CIP/SIP procedures between batches. This makes it ideal for biomanufacturers looking to avoid the time and costs associated with cleaning and validation.

The **MOTIV SU** system is particularly beneficial for processes involving chemicals that could contaminate future batches, as its replaceable, pre-built components ensure easy maintenance and quick readiness for reuse.

Built on the same platform as the **MOTIV**, the **MOTIV SU** offers a flexible, contamination-free solution that aligns with the growing trend of single-use systems in bioprocessing.

Now that you understand the different options for inline buffer formulation equipment, let's explore the specific features you need.






Features to **Consider**

As you shop for buffer formulation equipment, it's important to pay attention to the features because they directly impact the efficiency of your bioprocessing operations. **Here are essential features to consider:**



Automation



Flexibility



Usability



Accuracy and Precision



Blending Technology

Automation

Automation in buffer formulation systems allows for the consistent preparation of buffers across multiple batches, which is essential for maintaining product quality and reproducibility in bioproduction.

With automated systems, operators can trust the equipment to perform tasks precisely as programmed, minimizing the risk of human error, and ensuring that each batch meets strict specifications.

Incorporating 21 CFR Part 11 compliance into these automated buffer formulation systems further enhances the integrity and traceability of batch record keeping. By ensuring that all electronic records and signatures associated with buffer preparation are secure, validated, and auditable, manufacturers can confidently rely on their automated systems to produce consistent and high-quality buffers and maintain full regulatory compliance.

When every process step, from formulation to final product, is documented and verifiable, the manufacturer and end-user are protected from potential errors or discrepancies in production.

Flexibility

Flexibility allows a single buffer formulation system to efficiently produce large volumes of buffer, reducing the need for multiple systems and minimizing operational complexity. A flexible system also makes it easier to produce complex buffers that may require precise adjustments in pH, conductivity, or ionic strength.

Consider a flexible buffer formulation system that can adapt to different production needs, whether scaling up for large batches or fine-tuning conditions for specialized formulations. This adaptability ensures the equipment can meet the diverse and evolving demands of bioproduction, which ultimately maintains efficiency, reduces downtime, and simplifies the overall process of generating buffers

Usability

Usability is a key focus in the design of Asahi Kasei Bioprocess's buffer formulation equipment. Our "Built for You" approach emphasizes user functionality, safety, and comfort. By prioritizing ergonomic design, we ensure that our equipment is easy to operate and reduces user strain, enhancing productivity and efficiency on the floor.

A user-friendly interface is equally important in buffer formulation equipment because it ensures that operators can easily navigate and control the system. Operators may struggle to achieve the desired results if the automation software is complex or difficult to use



Asahi Kasei Bioprocess equips each MOTIV system with OCELOT™ System Control, a control software that integrates and interfaces with your plant-wide control system.

OCELOT offers intuitive interfaces with features like drag-and-drop functionality, making it accessible even to operators with less specialized training.

Accuracy & Precision

Accuracy and precision are largely driven by the high-quality sensors incorporated into buffer formulation systems. These sensors provide continuous, real-time feedback on the buffer's properties, such as pH and conductivity, ensuring that the buffer remains within the specified parameters.

With such precise control, the equipment can be programmed to maintain buffer conditions within strict limits. If the buffer's parameters deviate from the set range, the system can either shut down or quickly correct itself to prevent out-of-spec conditions. This high level of accuracy minimizes waste, as it ensures that the buffer remains usable and in specification for extended periods, whether minutes or hours, depending on the operational needs.

Blending Technology

Blending technology ensures buffers are mixed efficiently and precisely. One key aspect is recirculation blending technology, which continuously circulates the buffer to achieve a homogeneous mixture. This technology ensures that the blending process is thorough without being too aggressive, avoiding potential damage to sensitive materials.

The recirculation process allows the buffer to be blended to a rough state initially, after which sensors monitor the mixture to fine-tune its composition.

The blending technology and sensors adjust the mixing parameters to achieve the desired specifications. This real-time feedback loop ensures that the final buffer meets the required standards, enhancing accuracy and quality.

Speed is another important factor in blending technology. Efficient systems can blend the buffer and achieve the desired specifications rapidly, often within a minute. The recirculation blending technology facilitates this quick turnaround by continuously mixing the buffer and using pressure and sensor data to navigate the blending process effectively.

In addition to the essential buffer formulation equipment features, remember to consider your budget and the available vendor options. Let's explore these aspects in more detail.



Budget & Vendor Considerations

Before purchasing buffer formulation equipment, evaluate your budget and vendor options. Your budget helps ensure your investment is financially feasible and delivers optimal value. At the same time, selecting the right vendor affects the quality, reliability, and support you receive. Thoughtful consideration of these factors will position you for long-term success.

Budget

There are four core parts of your budget to consider when purchasing equipment for buffer formulation:



Initial investment

Operating costs

Maintenance and support

Trade-offs

Initial Investment

Your initial investment has a long-term impact on operational efficiency and cost-effectiveness. Even with a higher initial price, investing in higher-quality buffer formulation equipment can provide significant long-term benefits.

Premier buffer formulation equipment requires less maintenance and experiences reduced downtime, increasing productivity and cost savings over time. Since consistency and reliability are critical in buffer formulation, minimizing downtime with top-tier equipment is essential for smooth and efficient operations.

Operating Costs

Consumables such as **MOTIV® SU** tube sets help to maintain buffer purity and ensure the efficiency of each batch. These consumables must be replaced regularly to maintain optimal performance and avoid contamination, which contributes to the equipment's ongoing operational costs. By accounting for the cost of these consumables, you can more accurately estimate the total cost of ownership over the equipment's lifespan.



Preventative Maintenance & Support

Preventative maintenance is key to keeping your buffer formulation equipment running smoothly, helping to avoid unexpected failures and reduce downtime. Quick and efficient support from your vendor, whether remotely or on-site, is vital in resolving any issues and minimizing the impact on your operations.

Access to expert support is essential, as it gives operators the confidence to address technical challenges swiftly and effectively.

Consistent, reliable support ensures the longevity of your equipment and maintains high levels of productivity and cost-efficiency in your processes.

Trade-offs

While the initial investment in buffer formulation equipment may be substantial, it can lead to significant cost savings in other areas, particularly storage. Traditional buffer preparation methods require large volumes of pre-formulated buffers to be stored in massive tank farms, which not only demands substantial space but also incurs ongoing costs related to the maintenance, temperature control, and security of these storage facilities.

By using advanced buffer formulation equipment, you can produce buffers on demand and eliminate the need to store thousands of liters of prepared buffer solutions.

This frees up valuable space while reducing the financial burden associated with maintaining large storage infrastructures. Plus, on-demand buffer formulation minimizes the risk of buffer degradation over time, ensuring that the buffers used in your processes are always fresh and within the required specifications.

Vendors

When you're evaluating vendors, keep these considerations in mind:

- ▶ Customer support and training
- ▶ Reputation and experience
- ▶ Warranty and service agreements
- ▶ Regulatory compliance

Customer Support & Training

A vendor with dedicated customer support ensures researchers receive the help they need throughout the buffer formulation equipment's lifecycle. This begins with comprehensive training programs that enable operators to use and maintain the equipment effectively.

By offering hands-on training, educational materials, and recorded tutorials, the vendor makes it easier for users to implement and operate the equipment smoothly. Responsive customer support teams are also crucial, providing quick solutions to technical issues and queries.

On-site setup assistance and training that enhances the user experience will also instill confidence in operators and ensure successful equipment use.

Building correct buffer recipes is another critical part of optimizing your buffer formulation process. Expert guidance from your vendor allows you to navigate the intricacies of recipe development and ensure you understand the specific parameters necessary for your formulation.

Whether it's selecting the appropriate concentrations, adjusting pH levels, or achieving the perfect balance of components, knowledgeable support staff can provide the insights needed to refine your recipes.





Reputation & Experience

A strong reputation is founded on a history of reliability, innovation, and customer satisfaction. Selecting a reputable vendor significantly lowers the chances of facing challenges like product defects, delays, or insufficient support—factors that can disrupt research timelines and compromise results.

Industry recommendations and positive feedback from colleagues can offer valuable insights into a vendor's reliability, product quality, and customer service.

Therefore, it's wise to choose a vendor with a proven track record of consistently meeting customer expectations.

Warranty & Service Agreements

Warranty and service agreements are crucial when choosing a vendor for buffer formulation equipment. They offer essential support and assurance throughout the equipment's lifecycle. The warranty, typically lasting one to two years, protects against defects and provides confidence in the equipment's quality.

Service agreements can be customized to match the buyer's needs. Highly skilled teams might opt for basic coverage, while those needing more support may choose comprehensive agreements. These agreements ensure prompt technical assistance, troubleshooting, and preventative maintenance, which reduces downtime and maximizes equipment performance.

Regulatory Compliance

Choose a vendor that prioritizes regulatory compliance. Adherence to industry standards like Good Manufacturing Practices (GMP), Good Laboratory Practices (GLP), and 21 CFR Part 11 (for automation) ensures equipment meets strict quality and safety requirements, helping to reduce the risk of contamination or errors.

Vendors with ISO certification are committed to high-quality standards and effective quality control.

Validation and documentation support, such as that offered by Asahi Kasei Bioprocess, helps customers meet regulatory requirements and prepare for audits.

Certifications like UL, ASME, and ATEX further ensure equipment safety and suitability for use in hazardous environments, providing confidence in regulatory compliance and protecting the integrity of your processes.

In summary, investing in equipment that aligns with your budget and is sourced from a reputable vendor improves the efficiency and compliance of your buffer formulation process.

Conclusion

As illnesses become more complex and widespread, so does the need for effective biopharmaceuticals at scale. Large-scale bioproduction provides a path forward, and inline buffer formulation plays a key role in meeting the growing demand by enabling more flexible, efficient, and scalable operations.

Inline buffer formulation requires reliable, high-performing equipment to ensure consistent buffer preparation, minimize operational downtime, and maintain the stringent quality standards necessary for large-scale biopharmaceutical production.

When investing in inline buffer formulation equipment, consider the full budget, including initial investment, operating costs, maintenance, and trade-offs. While upfront costs may be high, reliable equipment offers long-term savings through increased efficiency, reduced downtime, and minimized storage needs.

Selecting the right vendor is crucial for success. A reputable vendor ensures reliability and innovation, while strong support and training help your team operate the equipment effectively. Comprehensive warranties and service agreements minimize downtime, and choosing a vendor with strict regulatory compliance ensures your equipment meets industry standards.

*Asahi Kasei Bioprocess offers innovative inline buffer formulation equipment, including the **MOTIV[®] 3**, **MOTIV[®] 5**, **MOTIV[®] plus**, and **MOTIV[®] SU**. Our products are built for you to ensure all your needs are met. To learn more about our products, [contact us](#).*



AsahiKASEI

BIOPROCESS

Contact Information

North and South America

Glenview, IL, USA
Tel: +1-847-556-9700
Fax: +1-847-556-9701

Email: info.us@ak-bio.com

Europe

Cologne, Germany
Tel: +49-221-995007-59
Fax: +49-221-9950077-59

Email: info.eu@ak-bio.com

India and Southeast Asia

Mumbai, India
Tel: +91-22-6710-3962
Fax: +91-22-6710-3979

Email: info.jp@ak-bio.com

East Asia and Oceania

Glenview, IL, USA
Tel: +1-847-556-9700
Fax: +1-847-556-9701

Email: info.jp@ak-bio.com

China

Shanghai, China
Tel: +86-(0)21-6391-6111
Fax: +86-(0)21-6391-6686

Email: info.jp@ak-bio.com