



Saving Costs Without Cutting Corners: Seal Replacement as a Cost Avoidance Strategy

Chuck Haeger
*Production And Order
Fulfillment Manager
Asahi Kasei Bioprocess*

With more than 35 years of experience in production management – including 15 with an emphasis on assembly/qualification of purification and synthesis columns – Chuck Haeger has a unique brand of expertise he leverages to empower his teams to deliver high quality products, all while also managing spare parts inventory and order fulfillment.

Rather than being a nonessential cost, routine seal replacement is a strategic cost avoidance tactic...

In biopharmaceutical manufacturing, cost-saving measures are often synonymous with risk—especially when they compromise product quality, process integrity, or regulatory compliance. However, not all savings come at a cost. Proactively replacing worn or degraded seals in synthesis columns is a small operational investment that can prevent costly failures, downtime, and even batch loss. Rather than being a nonessential cost, routine seal replacement is a strategic cost avoidance tactic—one that upholds quality, ensures process consistency, and safeguards long-term profitability. This article explores how a proactive seal management program can deliver real ROI, and why this seemingly minor maintenance activity is essential for operational excellence in regulated manufacturing environments.



COMPARING THE COST OF SEAL REPLACEMENT VS. COST OF CONTAMINATION

While replacing seals on synthesis columns may seem like a routine maintenance task, the cost differential between proactive replacement and the fallout from a contamination event is stark. Routine replacement of column seals is a low-cost, high-impact strategy that mitigates the risk of product contamination and regulatory non-compliance. Seals fabricated from materials such as PTFE, FKM, or EPDM often represent a considerable investment depending on chemical exposure and pressure requirements. Replacement during scheduled preventive maintenance adds minimal downtime and can be integrated into existing equipment validation protocols. In contrast, seal failure can lead to undetected ingress of contaminants or cross-contamination between synthesis steps—events that carry significant financial and regulatory consequences.



Saving Costs Without Cutting Corners: Seal Replacement as a Cost Avoidance Strategy


Chuck Haeger
*Production And Order
Fulfillment Manager
Asahi Kasei Bioprocess*

With more than 35 years of experience in production management – including 15 with an emphasis on assembly/qualification of purification and synthesis columns – Chuck Haeger has a unique brand of expertise he leverages to empower his teams to deliver high quality products, all while also managing spare parts inventory and order fulfillment.

By adhering to lifecycle-based replacement schedules and material compatibility assessments, facilities can prevent deviations, maintain product integrity, and ensure audit readiness.

In a 2022 FDA warning letter to a U.S. pharmaceutical manufacturer, investigators cited the failure to routinely inspect and replace degraded gaskets and seals in critical process equipment, resulting in visible black particulates in finished drug product (FDA Warning Letter 624650)¹. The affected batches were rejected, and the firm was required to halt production pending remediation and revalidation. In similar documented industry cases, the total cost of a single contamination event—including rejected product, investigation, batch recalls, and production downtime—can exceed hundreds of thousands to millions in damages.

From a risk-based manufacturing perspective, the ROI on proactive seal replacement is clear. By adhering to lifecycle-based replacement schedules and material compatibility assessments, facilities can prevent deviations, maintain product integrity, and ensure audit readiness. This is not just about maintenance—it's about protecting process reliability, patient safety, and the bottom line.



ROI OF PROACTIVE MAINTENANCE IN BATCH-BASED SYNTHESIS

In batch-based synthesis operations, proactive maintenance isn't just a best practice—it's a measurable contributor to return on investment (ROI). Unlike continuous processes where interruptions can be buffered by upstream redundancy, batch manufacturing is acutely vulnerable to equipment downtime and quality deviations. By contrast, implementing a structured preventive maintenance program—targeting critical components such as column seals, valves, and instrumentation—can dramatically reduce unplanned outages. For example, regular inspection and replacement of synthesis column seals at validated intervals may cost \$3,000–\$5,000 annually per line, including parts and labor. Yet, this investment can prevent contamination or equipment failures that could result in batch rejections, production stoppages, or extended deviation investigations—each with costs easily exceeding \$100,000.

¹ <https://www.fda.gov/inspections-compliance-enforcement-and-criminal-investigations/warning-letters/ultra-seal-corporation-624650-03142022>.



Saving Costs Without Cutting Corners: Seal Replacement as a Cost Avoidance Strategy

Chuck Haeger

Production And Order
Fulfillment Manager
Asahi Kasei Bioprocess

With more than 35 years of experience in production management – including 15 with an emphasis on assembly/qualification of purification and synthesis columns – Chuck Haeger has a unique brand of expertise he leverages to empower his teams to deliver high quality products, all while also managing spare parts inventory and order fulfillment.

For manufacturers looking to improve reliability while keeping costs in check, AI-powered seal monitoring offers a clear path to smarter, more efficient operations.

Moreover, facilities that proactively maintain equipment are less likely to incur regulatory citations, thereby reducing compliance risk and the potential costs of revalidation or remediation. When calculated over a fiscal year, the avoided costs from just two or three averted irregularities often exceed the full annual maintenance spend, delivering a strong positive ROI. Ultimately, in the high-stakes environment of Good Manufacturing Practices (GMP)-regulated batch production, proactive maintenance is not an overhead—it's a form of operational insurance with predictable returns.



AI-ENHANCED SEAL MANAGEMENT IN BIOMANUFACTURING: COST-EFFECTIVE, PREDICTIVE, AND PROACTIVE

Artificial Intelligence (AI) is rapidly gaining recognition as a key driver of autonomous manufacturing across industries. While it has yet to be fully integrated into seal maintenance practices within biomanufacturing, AI holds tremendous potential to revolutionize how manufacturers manage seals in synthesis columns—shifting maintenance from reactive, schedule-based routines to proactive, data-driven strategies. By analyzing historical process data—such as pressure patterns, temperature fluctuations, and past seal failure timelines—AI algorithms can identify subtle trends and anomalies that indicate early-stage seal wear. This enables facilities to predict when a seal is likely to fail, allowing for timely replacement that avoids both premature servicing and costly unplanned downtime.

In advanced implementations, AI systems can even correlate seal degradation patterns with batch performance, offering deeper insight into how seal health impacts yield and product quality. Integrating AI into maintenance workflows doesn't necessarily require massive infrastructure upgrades; many platforms can operate alongside existing sensors and control systems, making the transition both accessible and scalable. For manufacturers looking to improve reliability while keeping costs in check, AI-powered seal monitoring offers a clear path to smarter, more efficient operations.



Saving Costs Without Cutting Corners: Seal Replacement as a Cost Avoidance Strategy

Chuck Haeger

Production And Order
Fulfillment Manager
Asahi Kasei Bioprocess

With more than 35 years of experience in production management – including 15 with an emphasis on assembly/qualification of purification and synthesis columns – Chuck Haeger has a unique brand of expertise he leverages to empower his teams to deliver high quality products, all while also managing spare parts inventory and order fulfillment.

Recent research in bioprocessing demonstrates how machine learning can significantly improve maintenance strategies for critical components like seals in synthesis columns.

Recent research in bioprocessing demonstrates how machine learning can significantly improve maintenance strategies for critical components like seals in synthesis columns. A study applied advanced algorithms—Principal Component Analysis (PCA), Isolation Forests (IF), and Long Short-Term Memory (LSTM) autoencoders—to Protein A chromatography systems and successfully detected structural anomalies and integrity issues up to four cycles before failure.² Although the focus was column performance, many of these anomalies stem from seal deterioration or compromised packing, which manifest as subtle shifts in pressure, flow, or resin performance. By leveraging these early warning signals, manufacturers can transition from fixed-interval seal replacements to condition-guided maintenance, replacing seals only when needed. This reduces both material costs and labor while ensuring uninterrupted operation. In essence, AI-powered monitoring offers a low-investment, high-impact solution to bolster equipment reliability and drive down total cost of ownership in biomanufacturing.



In batch-based synthesis, the cost of replacing seals pales in comparison to the financial and reputational damage caused by contamination or process failure. What may seem like a minor maintenance task is, in fact, a pivotal safeguard—preserving product integrity, minimizing unplanned downtime, and ensuring regulatory compliance. By investing in routine seal replacement, manufacturers aren't simply maintaining equipment; they're making a strategic decision to avoid far greater costs down the line. In today's margin-conscious environment, this kind of proactive maintenance is not an expense—it's a smart, measurable investment in long-term operational resilience.

Questions?

Email me at
Chuck.Haeger@ak-bio.com



² <https://www.sciencedirect.com/science/article/abs/pii/S0021967322006781>.