

DAC Ergo Package

Safe and Efficient Maintenance of
Large-Scale DAC Columns (≥ 60 cm i.d.)



Features:

- » Three modes: process, travel, service
- » Quick-release bottom plate assembly
- » Transition from process mode to service mode in approximately 2 minutes

This document is confidential and solely owned by Asahi Kasei Bioprocess America, Inc. The use of this document, whole or part, or the facsimile of, is expressly forbidden without the written authorization of an AKBA officer.

+1.847.556.9700 | fluidmgmt.ak-bio.com

AsahiKASEI
BIOPROCESS