


MOTIV™ 5

Making the Complex Simple

Our 5-pump inline buffer formulation systems can produce more complex buffers constructed from up to five streams. Our systems provide the ultimate flexibility to generate mixtures using a combination of conductivity, pH and mass flow control. As with all of our inline buffer systems, OCELOT™ System Control will automate the process to ensure accuracy and repeatability.



powered by
OCELOT
SYSTEM CONTROL

MOTIV 52	MOTIV 512	MOTIV 520
20 – 200 L/h	120 – 1200 L/h	200 – 2000 L/h
Mobile console format	Space-saving skid format	
Process development	Clinical/large-scale manufacturing	
Up to 20X dilution	Up to 20X dilution	
770 mm W x 1360 mm D x 1600 mm H	1530 mm W x 1070 mm D x 1830 mm H	

Note: Dimensions are approximate.
 Customized systems for unique applications are also available.

[REQUEST A QUOTE >](#)